


ahn myTube MCT

PRODUCT DATA SHEET

Micro Tube 1.5 mL



Product Description

Catalogue Number	3-205-80-7
Product Description	AHN myTube MCT, 1.5 mL, black, bag
Sales Unit	Case / 16 x 500 tubes

Product Characteristics

Graduations	Yes
Sterilised	No
Autoclavable (Max °C, Period)	Yes, 121°C, 15 min
Centrifuge Specification	Up to 25.000 x g
Temperature Range	-80°C to +121°C
Purity	Free of RNase, DNase, Human DNA, Pyrogens, PCR Inhibitors

Size

Max. Volume	1.7 mL
Height of Product	40.9 mm

Material / Colours

Product Material	Polypropylene (PP)
Colour of Product	Clear

Packaging

Type of Packaging	Bag
Pieces / Bag	500 tubes
Pieces / Case	8000 tubes
Dimension Case Carton (H x W x D)	400 x 600 x 400 mm

Storage and Handling

Temperature	Room Temperature, ideally 24°C
UV Light	Protect from UV light

Note: Improper storage can initiate degradation, which results in colour changes and can have negative effects on the physical properties of the products.

Safety: This product is not classified as dangerous.

