

PRODUCT DATASHEET

AHN myTube® Centrifuge Tube CT



Product Description

Catalogue number	3-281-25-8
Product description	AHN myTube® CT 50 mL, clear/grey, self-standing, sterile, bag
Sales unit	Case / 20 x 25 tubes
GTIN sales unit	4255662801335

Product Characteristics

Graduations	Yes
Sterilised	Yes
Autoclavable (max. °C, period)	Yes, 121 °C, 15 min
Centrifuge specifications	Up to 7.000 x g
Temperature range	-80°C to +121°C
Purity	Free of RNase, DNase, Endotoxin, Pyrogens, non cytotoxic

Size

Max. Volume	50 mL
Height	115.0 mm

Material / Colours

Product material	Polypropylene (PP)
Colour of product	Clear/grey

Packaging

Type of packaging	Bag
Pieces / Bag	25 tubes
Pieces / Case	500 tubes
Dimension case carton (H x W x D)	355 x 600 x 400 mm
Weight	8,45 kg

Packaging

Temperature	Room temperature, ideally 24°C
UV light	Protect from UV light

Note	Improper storage can initiate degradation, which results in colour changes and can have negative effects on the physical properties of the products.
Safety	This product is classified as not dangerous.

Technical Drawing

